

ICPS-8100

Shimadzu
Sequential Plasma Spectrometer





**Now Offered with Ultimate Improvements to Satisfy
All Requirements for ICP Emission Analysis**

ICPS-8100

Shimadzu ICPS-8100 Twin Sequential High Frequency Plasma Emission Spectrometer SEQUENTIAL PLASMA SPECTROMETER


Since developing its first ICP emission spectrometer in 1977, Shimadzu has sold over 1200 multi-channel and sequential model ICP emission spectrometers and has remained a world leader in ICP technology. The ICPS-8100 is a top-of-the-line ICP emission spectrometer that not only incorporates the accumulated experience, performance, and quality of Shimadzu ICP emission spectrometers, but also features significantly improved performance, reliability, and ease-of-operation.

Advances in High Frequency Plasma

By introducing a carrier gas (Ar) and sample solution into plasma generated by passing a high-frequency current through an induction coil, a donut-shaped plasma with a lower temperature at the center than the perimeter is formed. This plasma is called inductively coupled plasma (ICP). Emission spectrometers that use this plasma as a light source offer many advantages, including a limitless wide range of applicability, over elemental analysis systems that use conventional light sources.

The ICPS-8100 minimizes losses by supplying a maximum 1.8 kW high-frequency output to the high-frequency coil via proprietary Shimadzu circuitry. Since a 27.12 MHz frequency is used to increase thermal energy, elemental excitation efficiency is high as well, resulting in an extremely high emission intensity. This allows introducing all types of solvents, such as organic solvents, hydrofluoric acid, or saturated salt water.

■ The twin sequential method provides stable analysis.



■ Benefits of the ICP light source, which is capable of highly accurate analysis of a wide variety of elements in a wide variety of samples in a wide range of concentrations

High detection capability to ppb-levels

Low interference between elements

Wide analytical concentration range

Outstanding stability and reproducibility

Contents

P 04 - Features

P 08 - Analytical Operations

P 15 - Specifications

P 06 - System Contents

P 12 - Sample Injection System

P 16 - Installation Requirements

P 07 - Application Fields

P 14 - Peripheral Equipment

Advanced Resolution, Speed, and Stability

Each of the two high-performance sequential monochromator units, capable of both high resolution and high speed, is equipped with a dedicated computer to achieve faster measurement and data processing speeds.

The monochromators use highly reliable Czerny-Turner mountings to provide sharp spectral lines and stability.

High Speed Measurement


Quantitative analysis and calculation of semi-quantitative values for **72 elements in 3 minutes**

Wide Wavelength Range

Vacuum monochromators and wide wavelength range from **160 to 850 nm**


High Resolution

Unrivalled **0.0045 nm** high resolution



Trace elements in pure tungsten (0.0045 nm resolution)

Measurement Solution:
1 g/100 mL + Fe, Co, and Cr in W



Wavelength profile of Zn 209.191 and As 189.042 in iron (0.0045 nm resolution)


Measurement Solution:
1 g/100 mL + Zn and As in high-purity Fe



Top Grade ICP


Comparison of Resolution 1 Sm in Ce

Measurement Solutions: Ce 1000 ppm and
Ce 1000 ppm + Sm 2 ppm



Comparison of Resolution 2 B in Steel

Measurement Solutions:
1 g/100 mL + B in high-purity Fe



Importance of Wavelength Resolution

High Resolution Required for Analysis of Metals, Rare Earths, and Rocks

High resolution is essential for analysis of metals, rare earths, and rocks. Improving the resolution helps improve separation between the target measurement wavelength and wavelengths of interfering elements or primary constituent elements. This enables highly accurate analysis of target elements down to trace levels, without being affected by interfering or primary constituent elements.




Twin Sequential Monochromators Provide Higher Analytical Efficiency and Ensure Stability for Long Periods

Twin spectrometer systems include two sequential monochromators that can be used for a wide variety of analytical applications. The two monochromators can be controlled independently for efficient analysis of samples. The software automatically selects the optimal monochromator for the wavelength being measured and sets the various analytical conditions. Monochromator 1 is always kept under vacuum conditions. Since the condensing lens can be cleaned while maintaining a vacuum, the

interior of the monochromator is not exposed to atmospheric conditions.

The vacuum monochromator enables high-precision analysis of elements that have highly sensitive analysis lines in the vacuum ultraviolet region, such as P, S, B, I, Br, and Al.


Since gas purging is not necessary, there is no contamination or variability caused by purge gas convection and the stabilization time is short. Therefore, it enables stable analysis over long periods.



Provides Highly Accurate Analytical Results for All Kinds of Fields, Such As the Detection of Ultra Trace Elements and Analysis of Chemical Compositions

The Shimadzu ICPS-8100 High Frequency Sequential Plasma Emission Spectrometer can be used wherever high accuracy and a broad scope of analytical assessment is required, whether analyzing ultra trace elements or high concentrations in chemical composition analysis. Applications include elemental analysis for research and

development, analysis of important elements for production control, and analysis for environmental management, such as water quality control, which has a significant impact on our lives.



ICPS-8100 Analytical Sensitivity Chart

Legend: Below 1 ppm 1 to 10 ppb 10 to 100 ppb Other

1a	2a	3b	4b	5b	6b	7b	8	1b	2b	3a	4a	5a	6a	7a	0		
H ¹															He ²		
Li ³	Be ⁴									B ⁵	C ⁶	N ⁷	O ⁸	F ⁹	Ne ¹⁰		
Na ¹¹	Mg ¹²									Al ¹³	Si ¹⁴	P ¹⁵	S ¹⁶	Cl ¹⁷	Ar ¹⁸		
K ¹⁹	Ca ²⁰	Sc ²¹	Ti ²²	V ²³	Cr ²⁴	Mn ²⁵	Fe ²⁶	Co ²⁷	Ni ²⁸	Cu ²⁹	Zn ³⁰	Ga ³¹	Ge ³²	As ³³	Se ³⁴	Br ³⁵	Kr ³⁶
Rb ³⁷	Sr ³⁸	Y ³⁹	Zr ⁴⁰	Nb ⁴¹	Mo ⁴²	Tc ⁴³	Ru ⁴⁴	Rh ⁴⁵	Pd ⁴⁶	Ag ⁴⁷	Cd ⁴⁸	In ⁴⁹	Sn ⁵⁰	Sb ⁵¹	Te ⁵²	I ⁵³	Xe ⁵⁴
Cs ⁵⁵	Ba ⁵⁶	L ⁵⁷	Hf ⁷²	Ta ⁷³	W ⁷⁴	Re ⁷⁵	Os ⁷⁶	Ir ⁷⁷	Pt ⁷⁸	Au ⁷⁹	Hg ⁸⁰	Tl ⁸¹	Pb ⁸²	Bi ⁸³	Po ⁸⁴	At ⁸⁵	Rn ⁸⁶
Fr ⁸⁷	Ra ⁸⁸	A ⁸⁹															
L	La ⁵⁷	Ce ⁵⁸	Pr ⁵⁹	Nd ⁶⁰	Pm ⁶¹	Sm ⁶²	Eu ⁶³	Gd ⁶⁴	Tb ⁶⁵	Dy ⁶⁶	Ho ⁶⁷	Er ⁶⁸	Tm ⁶⁹	Yb ⁷⁰	Lu ⁷¹		
A	Ac ⁸⁹	Th ⁹⁰	Pa ⁹¹	U ⁹²	Np ⁹³	Pu ⁹⁴	Am ⁹⁵	Cm ⁹⁶	Bk ⁹⁷	Cf ⁹⁸	Es ⁹⁹	Fm ¹⁰⁰	Md ¹⁰¹	No ¹⁰²	Lr ¹⁰³		

Detection Limit of ICPS-8100 (ppb)

- Below 1 ppm
- 1 to 10 ppb
- 10 to 100 ppb
- Other




Analytical Operations Performed and Instrument Controlled Easily by a PC

All analytical and instrument control operations are performed by a PC. Parameters for instrument control and analysis can also be easily changed.

Each monochromator includes a dedicated control computer that allows confirming analytical conditions, checking calibration curves, or processing data and profiles from samples already measured, even while analyses are in progress. Respective processes can be performed without interrupting analyses, which can significantly shorten analysis times.

All operations can be performed via software windows, from igniting the torch and confirming the plasma status, to confirming the safety status or controlling the instrument.

Instrument Start up



Top Grade ICP

Qualitative Analysis of All Elements

When the group is selected, the measurement procedure is displayed.



Elements are selected in advance.

High Accuracy Qualitative Analysis

When the group is selected, the measurement procedure is displayed.

Simply select the elements to be measured.

No.	Type	Repeat	Sample Name	Weight	Table	Status
1	Sample	1	BLANK	0.0000	1	
2	Sample	1	Sample 1	0.0000	2	
3	Sample	1	Sample 2	0.0000	3	
4	Sample	1	Sample 3	0.0000	4	
5	Sample	1	STD	0.0000	5	
6						
7						
8						



Item	Value	Range	Step
Zn	12		
B	1.9		
Al	1.1		
Fe	5.3		
Al	20		
Ge	30		
Hg	13		
Li	0.002		
Fe	0.00		
Pd	0.63		

Item	Value	Range	Step
RF Power	1.2	0 - 1.0	
Observation Height	11.0	mm 11 - 25	2
Coolant Gas	14.0	L/min 6 - 20	2
Plasma Gas	1.20	L/min 0.4 - 1.4	0.2
Carrier Gas	0.70	L/min 0 - 1.55	0.05
Purge Gas	<input checked="" type="checkbox"/> ON		
Rinse Time			
Solvent Rinse-H	5	Solvent Rinse-L	15
Sample Rinse-H	5	Sample Rinse-L	45
Source	<input checked="" type="checkbox"/> Plasma	Chamber	<input checked="" type="checkbox"/> CH1
	<input type="checkbox"/> Hg-Lamp		<input type="checkbox"/> CH2
Eco. Run	<input type="checkbox"/> No		<input type="checkbox"/> Yes

Item	Value	Range	Step
Na	8.6		
Ni	5.4		
Zn	7.5		
B	1.9		
Si	1.4		
Co	10		
Cu	89		
Mg	0.03		
V	0.00		
Cr	0.46		
Mn	0.024		

Optimal measurement parameters are selected. Parameters can be set manually as well.

Automatically Sets Measurement Parameters for Qualitative to Quantitative Analysis

Quantitative Analysis




No.	Element	Type	Number	Wavelength	Integ line	ATT	JAMP	Condition	Delta	Menu
1	B	Sequential	5	243.772	50	7	x10	1	2M	
2	Na	Sequential	11	588.995	50	9	x10	1	25	
3	Al	Sequential	13	396.153	50	8	x10	1	2M	
4	Si	Sequential	14	251.612	50	6	x10	1	2M	
5	P	Sequential	15	178.267	50	9	x10	1	1M	
6	S	Sequential	16	180.731	50	9	x10	1	1M	
7	K	Sequential	19	766.491	50	15	x10	1	25	
8	Ni	Sequential	28	231.604	50	4	x10	1	1M	
9	Cu	Sequential	29	327.396	50	6	x10	1	2M	
10	Zn	Sequential	30	213.808	50	6	x10	1	2M	
11	As	Sequential								
12	Se	Sequential								

Integration time can be varied for each element.
All measurement parameters are set automatically.

Element	Sample	No.	Element	Conc.	Unit	Intensity	Type
STD1							
STD2		1	B	500000	ppm	.000000	
STD3		2	Na	500000	ppm	.000000	
STD4		3	Al	500000	ppm	.000000	

Element	Seq. Internal	Internal	Variable Internal	Fixed Internal
B (2M)			B (2M)	
Na (25)			Na (25)	
Al (2M)			Al (2M)	
Si (2M)			Si (2M)	
P (1M)			P (1M)	
S (1M)			S (1M)	

Type	Repeat	Sample Name	Weight	Table	Status
Standard	3	STD1	.000000	1	
Standard	3	STD2	.000000	2	
Standard	3	STD3	.000000	3	
Standard	3	STD4	.000000	4	
Sample	3	Sample1	.000000	5	
Sample	3	Sample2	.000000	7	
Sample	3	Sample3	.000000	8	




Capable of unmanned operation, including automatic measurement and automatic extinguishing of the torch.

Multiple calibration curves and enlarged views can be displayed with a single click.


Top Grade ICP

If qualitative analysis is performed, measurement parameters are automatically transferred for quantitative analysis. Analyses can be performed using optimal parameters from qualitative analysis results. Even if starting with quantitative analysis, **measurement parameters are set automatically.**

In addition, searching for interfering elements after measurements can be performed by **single-click operations on the enlarged peak view window.** The instrument automatically determines the optimal background position as well.



Search for and identify interfering peaks easily in the enlarged peak view



Enlarged views of profiles and searches for interfering elements can be displayed with single clicks.

Calibration method result - Manual				
Sample 1		Prev. Sample	Next Sample	
Sample Name	Sample 1			
Analysis date	2009/07/01 11:06			
<Intensity>				
Element	Mg	Ca	Fe	Zn
No. 1	539.964	4197.65	4.63065	1.20923
No. 2	541.555	4226.84	4.56807	1.21678
No. 3	540.164	4277.76	4.60248	1.21289
Average	540.561	4234.08	4.60040	1.21297
R	1.59052	80.1025	.062581	.007545
S	.866421	40.5395	.031342	.003773
CV	.160282	.957457	.681298	.311043
<Conc.>				
Element	Mg	Ca	Fe	Zn
No. 1	2.39016	12.3731	.359349	.073727
No. 2	2.39717	12.4595	.354103	.074258
No. 3	2.39104	12.6102	.356988	.073984
Average	2.39279	12.4809	.356813	.073990
R	.007017	.237065	.005246	.000531
S	.003822	.119978	.002627	.000266
CV	.159742	.961288	.736322	.358941










Measurement data is saved as a file. Of course, profile data is similarly saved. Data can be transferred to commercial software to freely create reports. Since the instrument is connected to the data processor via a LAN, data can be processed from a separate office.






















Profile data and measurement results can be pasted in commercial spreadsheet software to allow easy formatting of reports.

ICPS-8100 Accessories

Recommended Combinations for Sample Injection Systems

(Note) Purchase items individually for sample types without a part number (P/N)

Item	Nebulizers	Sample Take-Up Tubing
Standard Set ICPS-8100 Standard Accessories	Nebulizer, 10ES (P/N 046-00092-02) Designed for high-concentration samples, with high-efficiency nebulization. (Sample take-up rate approx. 1 mL/min) <ul style="list-style-type: none"> • Sample take-up tube ASSY, S-075 • Connector, QSM (P/N 046-00092-09) • Tube adaptor, 0735 (P/N 046-00092-10) These are included. 	Sample take-up tube ASSY, S-075 (P/N 046-00092-06) For "10ES" and "07ES" nebulizers. 
For Small Volume Samples	Nebulizer, 07ES (P/N 046-00092-01) A nebulizer with high-efficiency nebulization. (Sample take-up rate approx. 0.6 mL/min) <ul style="list-style-type: none"> • Sample take-up tube ASSY, S-075 • Connector, QSM (P/N 046-00092-09) • Tube adaptor, 0735 (P/N 046-00092-10) These are included. 	Sample take-up tube ASSY, S-075 (P/N 046-00092-06) For "10ES" and "07ES" nebulizers. 
High Salt Samples Used to inject high salt concentration samples	Nebulizer, 10ES (P/N 046-00092-02) (Sample take-up rate approx. 1 mL/min) <ul style="list-style-type: none"> • Sample take-up tube ASSY, S-075 • Connector, QSM (P/N 046-00092-09) • Tube adaptor, 0735 (P/N 046-00092-10) These are included. 	Sample take-up tube ASSY, S-075 (P/N 046-00092-06) For "10ES" and "07ES" nebulizers. 
Organic Solvent Samples Injection systems for organic solvents	Nebulizer, 10ES (P/N 046-00092-02) (Sample take-up rate approx. 1 mL/min) <ul style="list-style-type: none"> • Sample take-up tube ASSY, S-075 • Connector, QSM (P/N 046-00092-09) • Tube adaptor, 0735 (P/N 046-00092-10) These are included. 	Sample take-up tube ASSY, S-075 (P/N 046-00092-06) For "10ES" and "07ES" nebulizers. 
Water-Cooled Chamber Kit for Organic Solvent Samples Controls evaporation of solvents	<hr/>	<hr/>
Hydrofluoric Acid Sample Injection System (P/N 211-42853-03) Used for hydrofluoric acids	Nebulizer, 10CPS (P/N 046-00092-14) <ul style="list-style-type: none"> • PTFE Tube, 1.27×1500L (P/N 046-00092-03) • Tube adaptor, 0735 (P/N 046-00092-10) These are included. 	<hr/>

Chambers	Torches	Drains	Other
<p>Cyclone chamber, HE (P/N 046-00093-02)</p>  <ul style="list-style-type: none"> • Locking screw, 0152 (P/N 046-00093-92) • Seal, 0237 (P/N 046-00093-93) <p>These are included.</p>	<p>Torch (P/N 204-70272)</p> 	<p>Drain trap,8214 (P/N 046-00093-01)</p> <p>Drain trap for cyclone chamber (Indicated by "•")</p> 	<p>Ball joint clip (P/N 210-15508-01)</p> <p>Connects the cyclone chamber and torch.</p> <p>Clip (P/N 046-00993-01)</p> <p>Fixes the drain trap to the cyclone chamber.</p> <p>Drain tube (P/N 200-30864-24)</p> <p>This is the tube from the drain trap to the drain.</p>
<p>Cyclone chamber, HE (P/N 046-00093-02)</p>  <ul style="list-style-type: none"> • Locking screw, 0152 (P/N 046-00093-92) • Seal, 0237 (P/N 046-00093-93) <p>These are included.</p>	<p>Torch (P/N 204-70272)</p> 	<p>Drain trap,8214 (P/N 046-00093-01)</p> <p>Drain trap for cyclone chamber (Indicated by "•")</p> 	<p>Ball joint clip (P/N 210-15508-01)</p> <p>Connects the cyclone chamber and torch.</p> <p>Clip (P/N 046-00993-01)</p> <p>Fixes the drain trap to the cyclone chamber.</p> <p>Drain tube (P/N 200-30864-24)</p> <p>This is the tube from the drain trap to the drain.</p>
<p>Cyclone chamber, HE (P/N 046-00093-02)</p>  <ul style="list-style-type: none"> • Locking screw, 0152 (P/N 046-00093-92) • Seal, 0237 (P/N 046-00093-93) <p>These are included.</p>	<p>Torch for high-concentration salt solution samples (P/N 204-74323)</p> <p>Used when introducing an undiluted solution of a sample where the base element is of the % order of concentration.</p> 	<p>Drain trap,8214 (P/N 046-00093-01)</p> <p>Drain trap for cyclone chamber (Indicated by "•")</p> 	<p>Water bubbler (P/N 204-19281)</p> <p>Used to prevent blockages in the nebulizer when analyzing samples containing large amounts of sodium. Inserted in the carrier gas pathway.</p> 
<p>Chamber drain straight (P/N 211-80437-01)</p> 	<p>Organic solvent torch (P/N 204-77296)</p> <p>Used when analyzing organic solvent samples that are difficult to inject into the plasma.</p> 	<p>Drain kit (P/N 211-86140-91)</p> 	<p>Nebulizer holder kit type 2 (P/N 211-48062)</p> <p>Includes organic solvent resistant O-rings. They are attached to each chamber other than the cyclone chamber, and are used to support nebulizers.</p> 
<p>Water-cooled chamber kit (P/N 211-43472)</p> <p>This increases injection efficiency of organic solvent samples by cooling the chamber and suppressing sample evaporation within the chamber. *Cooling water circulator is necessary.</p> 	<p>Demountable torch (P/N 205-09627-01)</p> 	<p>Drain kit (P/N 211-86140-91)</p> 	<p>NCB-1200 (SP) (P/N 044-01910-01)</p> 
<p>Chamber ASSY (P/N 205-07778-02)</p> 	<p>Demountable torch (P/N 205-09627-01)</p> 	<p>Drain kit (P/N 211-86140-91)</p> <p>(Included in the hydrofluoric acid sample injection system)</p> 	

A Wide Variety of Peripherals

Accessories that enable automatic analysis and sample introduction at your discretion.

Autosampler AS-8T

(P/N 205-04940-02)

This type handles both beakers and test tubes. It can be loaded with 100 of the 20 mL test tubes or 50 of the 30 mL beakers. A table is necessary.

*Requires interface kit with table (P/N 211-80655).

Table size: W600xD600xH690 mm



AS-8T for organic solvents

(P/N 211-48059-01)

Hydrofluoric acid sample injection system HFS-2

(P/N 211-42853-03)


Samples composed primarily of silicates, such as rocks, soil, cement and ceramics, are insoluble or dissolve poorly in strong acids. When hydrofluoric acid is used to dissolve these samples, normal injection systems cannot be used as they are made of glass. A sample dissolved in hydrofluoric acid can be injected into the plasma directly by using an injection system made of fluorinated ethylene resin.



Peristaltic pump

(P/N 204-77310-02)

Used in the analysis of high viscosity samples. A fixed sample injection is possible.



Monochromator for internal standard method

(P/N 205-02165-02)

This is installed in the Shimadzu sequential plasma spectrometer to perform simultaneous internal standard analysis.

*Requires attachment kit (P/N 211-82923).

Specifications

Optical mount: 0.5 m Pachen-Runge mounting
 No. of grooves in the diffraction grating: 2700 grooves/mm
 Internal standard element: Y (371.0 nm)
 Light guide: by optical fiber
 Place of installation: Inside thermostatic chamber in Monochromator 1

Axial attachment AX-3

(P/N 211-41992)

*A cooling water system is necessary as a utility.

If water for the cooling system is not supplied from the public water system, a CA-1112 cooling water circulator and piping kit D are required.

Ultrasonic nebulizer UAG-1

(P/N 205-09295)

This is a sample injection system developed for high-sensitivity ICP analysis. Unlike the usual nebulizer that uses the negative pressure of the carrier gas, this nebulizer uses ultrasonic energy to nebulize a solution. Ultrasonic energy can produce large amounts of extremely fine particles. This enables the UAG-1 to perform analyses at a high sensitivity of 10 to 100 times the sensitivity of conventional nebulizers.

Size: Main unit W320xD380xH540 mm
 Power supply: W382xD360xH159 mm
 Power supply: Single phase 100 V, 50/60 Hz, 10 A

*A cooling water system is necessary as a utility.

If water for the cooling system is not supplied from the public water system, a CA-1112 cooling water circulator and piping kit C are required.



Table ND with caster

(P/N 219-96005) (included) Size: W600xD480xH750 mm

Hydride generator HVG-ICP

(P/N 211-40981)

The elements within the sample are reduced and vaporized by the nascent hydrogen generated in the decomposition of sodium borohydride. Only the gas phase is injected into the plasma to achieve measurement with a high degree of sensitivity. As, Se, Sn, Te, Bi etc. can be measured.

Size: W333xD210xH195 mm
 Power supply: Single phase 100 V, 50/60 Hz, 2 A
 *A table ND (with casters) is necessary.



Table ND with caster

(P/N 219-96005) Size: W600xD480xH750 mm

Cooling water circulator CA-1114

(P/N 044-01809-07)

Size: W354xD384xH851 mm 41 kg
 Power supply: Single phase 100 V, 50/60 Hz, 15 A

Piping kit C: UAG-1 only

(P/N 211-83633-03)

Piping kit D: AX-3 only

(P/N 211-83633-04)

Piping kit U: AG-1 and AX-3

(P/N 211-83633-05)



Low-temperature thermostatic water heater NCB-1200 (SP)

(P/N 044-01910-01)

For the water cooled chamber kit

Size: W210xD430xH639 mm 29 kg
 Power supply: Single phase 100 V, 50/60 Hz, 9.5 A



Specifications

1. Monochromator unit	
Monochromator 1	
Optical system	1 m Czerny-Turner mounting
No. of diffraction grating grooves	4960 grooves/mm
Wavelength range	160 to 372 nm
Reciprocal dispersion	0.15 nm/mm
Detector	Photomultiplier tube
Slit	Entrance slit 20 μm Exit slit 30 μm
Exhaust system	Rotary pump With oil backflow prevention valve
Temperature control	Available (± 0.1 °C)
Monochromator 2	
Optical system	1 m Czerny-Turner mounting
No. of diffraction grating grooves	4320 grooves/mm 1800 grooves/mm
Wavelength range	250 to 426 nm (4320 grooves/mm) 426 to 850 nm (1800 grooves/mm)
Reciprocal dispersion	0.17 nm/mm (4320 grooves/mm) 0.44 nm/mm (1800 grooves/mm)
Detector	Photomultiplier tube
Slit	Entrance slit 20 μm Exit slit 30 μm
Light guide	Optical fiber
Temperature control	Available (± 0.1 °C)
Spectrometer for internal standard method (optional)	
Optical system	0.5 m Pachen-Runge mounting
No. of diffraction grating grooves	2700 grooves/mm
Internal standard element	Y (371.0 nm)
Light guide	Optical fiber
Place of installation	Inside thermostatic chamber in Monochromator 1
Monochromator controller (installed on respective monochromators 1 and 2)	
CPU	With 32-bit microprocessor
Photometric unit (installed on each monochromator)	
Negative high-voltage power supply	Variable 16 steps
Photometric method	Sequential element measurement method
2. ICP light source	
Plasma torch	Quartz plasma torch
Nebulizer	Coaxial type, Pyrex glass
Spray chamber	Cyclone chamber, Pyrex glass
Torch vertical actuation mechanism	Computer controlled (11 mm to 25 mm)
3. Radio frequency generator	
Oscillator	Crystal oscillator
Frequency	27.120 MHz ± 0.05 % (ISM band)
Output	0.8, 1.0, 1.2, 1.4, 1.6, 1.8 kW
Output stability	Within ± 0.3 %
Radio frequency circuit element	Transistor
Control system	CPU control
Ignition method	Fully automatic ignition
With safety function	Radio frequency generator temperature fault detection

4. Four flow line gas controller (for plasma)	
Plasma gas	0 to 20.2 L/min steps
Auxiliary Gas	0 to 1.4 0.2 L/min steps
Carrier gas	0 to 1.55 0.05 L/min steps
Purge gas (Monochromator 1)	3.5 L/min
5. Data processor (software)	
OS	Windows XP/Windows Vista
No. of measurement wavelengths	Qualitative analysis 1: 72 elements \times 1 wavelength Qualitative analysis 2: 72 elements, maximum 216 wavelengths Quantitative analysis: 72 elements, maximum 72 wavelengths
Database	Analyzed wavelengths data: 72 elements, maximum 16 wavelengths Wavelength table: Approx. 110,000 wavelengths recorded
Analysis cards	100 cards
Quantitative analysis	Calibration curve sample, maximum 16 samples per card Drift correction Internal standard correction Background correction Blank signal elimination Matrix correction
6. Safety functions	
Plasma extinguishment detection	Stand temperature check
Argon gas pressure check	High-frequency power supply temperature check
Cooling water circulation check	Monochromator vacuum check
Stand door open/closed check	Monochromator temperature check
7. Autosampler AS-9	
Any sample is accessible by CPU control	
No. of samples	50
Sample container	20 mL test tube
Actuation method	X-Y actuation

8. External dimensions	
Units: mm	
W1990xD1040xH960 / Weight: Approx. 640 kg	

* Windows is a registered trademark in the United States and other countries of Microsoft Corporation.

* Intel and Celeron are registered trademarks of Intel Corporation in the United States and other countries or trademarks of the subsidiary company.


* Additionally, the company name and product name that has been described are the trademarks or registered trademarks of each company. * Neither TM nor [®] mark are described clearly in the text.

ICPS-8100 Installation Requirements

Prepare the following for installation.

Refer to the installation guidelines for details.

1. Installation room environment	Temperature: 18 to 28 °C
	Rate of temperature change: 2 °C/hr max.
	Humidity: 70 % max.
	Avoid places with a lot of vibration or dust.
	Heat generated from instrument is about 2200 kcal/hr (during 1.2 kW high frequency output, excluding heat discharged to exhaust duct)
2. Power source	3-phase: 200/220 V±10 %, 50/60 Hz, 20 A
	Single phase: 100 V±10 %, 50/60 Hz, 5 A
3. Grounding	Resistance: independent grounding below 30 Ω
4. Gas installation	Type: Argon gas of 99.95% or greater purity
	Adjust argon gas supply within 550 to 750 kPa pressure (maximum variation 100 kPa). Approximately one 7 m ³ gas cylinder is necessary for 5 hours of operation.
5. Cooling water	Prepare 12 L of distilled water as coolant for cooling high-frequency induction coil in main unit.
	If a UAG-1 ultrasonic nebulizer (optional) is included, one water cooling system is required.
6. Exhaust duct	Plasma stand / Exhaust gas is mostly argon, but also includes some metal vapors and solvent. Therefore, install exhaust fan blades for explosion resistant windows and exhaust ducts.
	High-frequency power supply / Ducting is not required, due to low heat generation.
7. License	Usage of this device needs to comply with the radio laws. A license for an installation using radio frequency needs to be obtained from the authority.



JQA-0376

- * Windows is a trademark of Microsoft Corporation of the United States of America.
- * TIFF is a trademark of Adobe Systems, Inc. of the United States of America.
- * All other company and product names indicated in this document are trademarks or registered trademarks of their respective companies.

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan
Phone: 81(3)3219-5641 Fax: 81(3)3219-5710

URL <http://www.shimadzu.com>

The contents of this brochure are subject to change without notice.