



**AAS SERIES**

**AA-7003 Automatic Flame/Graphite Furnace Atomic Absorption Spectrometer**



**EWAI** | SIMPLY  
PRECISE

**AA-7003 Automatic  
Flame/Graphite Furnace Atomic  
Absorption Spectrometer**



## New Features

- Excellent performance / price ratio
- Advanced built-in graphite furnace power supply
- Can be equipped with a graphite furnace autosampler
- Durable graphite tube, can be used for a long time at under 2000 °C
- Maximum 20 segment graphite furnace temperature program
- Powerful software that features easy-to-use method parameter setup and advanced calibration curve fitting
- Unique background reduction design that allows for accurate analysis of high salinity, high background samples such as soy sauce, food, etc.