



**AAS SERIES**

**AA-7050 Atomic Absorption Spectrometer**



# **AA-7050 Atomic Absorption Spectrometer**



## New Features

- The graphite furnace camera helps users locate injection probe, view sample status in graphite tubes and optimize furnace temperature programming.
- Graphite furnace gas saving mode prolongs the usage time of protection gas and therefore reduces the operating cost of the instrument.
- Graphite furnace double inner gas path permits an auxiliary gas in addition to the protection gas so as to improve the ashing process and obtain better results
- 6 coded lamp turret can automatically identify lamp element and record element lamp usage times.
- Deuterium lamp fine adjustment mechanism ensures perfect overlap of deuterium lamp and element lamp positions.
- Multi-element analysis function and AS-600 flame/graphite furnace auto sampler are equipped to realize fully automated analysis of multiple elements