



AAS SERIES
AA-7090 Atomic Absorption Spectrometer



AA-7090 Atomic Absorption Spectrometer



New Features

- AA-7090 flame and Zeeman graphite furnace tandem analyzer, with two background correction methods: longitudinal Zeeman effect correction and deuterium lamp correction
- Transverse heating, longitudinal Zeeman effect graphite furnace, introduced from Australia.
- Magnetic field can be varied from a low of 0.6 to a high of 1.1 Tesla in increments of 0.1 Tesla
- 8-Lamp automated lamp turret
- ESV
- Coded lamp recognition
- High efficiency super lamp power supply
- The burner head height can be automatically adjusted
- Programmable gasbox
- Protection gas saving system, high efficiency & reduce your cost