


## FTIR bench top spectrometer Interspec 200-X



### Specifications:

Wavelength range, IR	7000 to 400 $\text{cm}^{-1}$
Wavelength range, NIR	14000-5800 $\text{cm}^{-1}$ or 9000 -3850 $\text{cm}^{-1}$
Resolution, standard	1 $\text{cm}^{-1}$
Resolution, option	0.5 $\text{cm}^{-1}$
Beam at sample	10 mm dia.
Operating system	Windows based
Data acquisition system	18 bit, high speed
Interface	USB 2.0
Power	12 VDC, 30 W
Sample compartment	W20xD26xH16 cm
Dimensions	W59xD39xH19 cm
Weight	24 kg

**Interspectrum OU** - Toravere, 61602 Tartumaa, Estonia  
Tel.: +372 7383008, Fax: +372 7383065  
Email: [info@interspectrum.ee](mailto:info@interspectrum.ee)  
[www.interspectrum.ee](http://www.interspectrum.ee)