

Safety Covers

Preventing door glass explosion for user safety with a netted steel frame and polycarbonate sheet covering.

Models	Description	Cat. No.
OV-11	Safety Cover	AAA13601
OV-12	Safety Cover	AAA13602

Gaskets

Available very persistent viton (Fluorine rubber) door gasket upon request.

Models	Description	Cat. No.
OV-11	Silicone Gasket	FAA5146
	Viton (Fluorine rubber) Gasket	FAA5187
OV-12	Silicone Gasket	FAA5147
	Viton (Fluorine rubber) Gasket	FAA5188

Vacuum Pumps


Direct connecting oil sealed and backward flow preventive equipped rotary vacuum pumps, which are essentially used with vacuum ovens to get high-degree vacuuming condition. Oil mist trap to prevent oil inflowing is optional.

Models	Electrical Requirements	Cat. No.
TRP-6	230V, 50/60Hz	BEH73005K
	120V, 60Hz	BEH73006U
Related accessory	Oil Mist Trap	BEH73501

Cold Trap Baths

Preventing vacuum pumps in a condition of vacuuming as a cold trap to exclude water (vapor). Included 2 U-type glass traps.

Models	Electrical Requirements	Cat. No.
CTB-10	230V, 60Hz	AAH62011K
	230V, 50Hz	AAH62012K
Related accessory	U-type Glass Trap	AAA62501



Viton Gaskets



Vacuum Pumps

Mehrgo
Parto
Shar



Cold Trap Baths



Safety Covers

Mehrgo
Parto
Shar