

Mechanical Ultrasonic Cleaner

The Mechanical Series is our best-selling model. This low cost option offers powerful and reliable cleaning at an affordable price. All Qsonica cleaners are made with a stainless steel tank and chemically resistant plastic body. A cover is included with every unit.

A mesh basket or perforated tray must be ordered separately. A Basket is recommended for cleaning small items and the Tray is for use with heavier items or items that must be held upright and level.

Mechanical Ultrasonic Features

- 60-minute timer
- Continuous operation option
- Low cost and reliable
- 40kHz frequency

Mechanical Ultrasonic Specifications

Model #	Volume	Overall Size L x W x D	Tank Size L x W x D	Weight	Frequency
C50T	½ gallon (1.9L)	10 x 12 x 11.5" (254 x 305 x 292mm)	6 x 5.5 x 4" (152 x 140 x 102mm)	7 lbs.	40 kHz
C75T	¾ gallon (2.8L)	13.5 x 12 x 11.5" (343 x 305 x 292mm)	9.5 x 5.5 x 4" (241 x 140 x 102mm)	9 lbs.	40 kHz
C150T	1½ gallon (5.7L)	16 x 12 x 14.5" (406 x 305 x 368mm)	11.5 x 6 x 6" (292 x 152 x 152mm)	12 lbs.	40 kHz
C250T	2½ gallon (9.5L)	16 x 15.5 x 14.5" (406 x 394 x 368mm)	11.5 x 9.5 x 6" (292 x 241 x 152mm)	14 lbs.	40 kHz