

Power Adjustable Ultrasonic Cleaner

Power Adjustable Ultrasonic Cleaner Features

1. Display the parameters setting menu
2. Adjustable power :40%-100%
3. Digital setting ultrasonic cleaning time
4. Working time countdown display
5. Display actual temperature of the cleaning tank
6. Temperature setting Range: room temperature-80°C
7. Timing Range: 1-999 minutes adjustable
8. Equipped with a special stainless steel basket, reduce sound cover
9. The shell and cover are made by stainless steel
10. The equipment controlled by Single-chip controller
11. The latest working parameters can be stored in memory after power-off

1Page |



Power Adjustable Ultrasonic Cleaner Specifications

Model	Inner Tank	Volume	Frequency	power	Heating power	Drainage	Adjustable Temp.	Adjustable time	Rack
	L/W/H (mm)	L	KHZ	W	W		°C	min	
SB-120DTP	240×140×100	3	40	120	120	/	Ambient-80	1-999	√
SB-3200DTDN	300×155×150	6	40	180	120	/	Ambient-80	1-999	√
SB-5200DTDN	300×240×150	10	40	200	320	√	Ambient-80	1-999	√
SB-5200DTDN	300×240×150	10	40	250	320	√	Ambient-80	1-999	√
SB-5200DTDN	300×240×150	10	40	300	320	√	Ambient-80	1-999	√
SB-4200DTDN	300×240×200	14.4	40	400	320	√	Ambient-80	1-999	√
SB25-12DTDN	500×300×150	22.5	40	500	600	√	Ambient-80	1-999	√
SB25-12DTDN	500×300×150	22.5	40	600	600	√	Ambient-80	1-999	