

Ultrasonic Homogenizer Q55

The Q55 is a compact and cost effective ultrasonic processor that will occupy less bench space than any unit on the market.

This model is effective for standard cell disruption and many other small volume applications. Probes are available in three different sizes.



Ultrasonic Homogenizer Q55 Features

- Smallest unit available
- Thumb-switch or continuous operation
- Simple and effective operation

Ultrasonic Homogenizer Q55 Parts

- Generator
- Converter
- 1/8" diameter probe
- Power cable
- Converter cable
- Wrench set

Ultrasonic Homogenizer Q55 Features

Power Rating	55 Watts
Frequency	20 kHz
Dimensions	8" W x 7.5" L x 5.75" H
Voltage	110V, 50/60Hz



MPS
Mehrgan
Parto
Shar