

Wire Shelves

Stainless steel and electro polished shelf system is removed easily without using tools for easy cleaning.



Wire Shelves


Product group	Models	Description			Cat. No.
		Included	Maximum	W×L (mm/inch)	
Ovens	OF-02G	2	8	366×316 / 14.4×12.4	EDA8136
	OF-12G	2	12	446×366 / 17.6×14.4	EDA8135
	OF-22G	2	14	516×416 / 20.3×16.4	EDA8134
	ON-02G	2	8	366×316 / 14.4×12.4	EDA8136
	ON-12G	2	12	446×366 / 17.6×14.4	EDA8135
	ON-22G	2	12	576×419 / 22.7×16.5	EDA8137
	OF/ ON-01E	2	10	366×316 / 14.4×12.4	EDA8136
	OF/ ON-11E	2	14	446×366 / 17.6×14.4	EDA8135
	OF/ ON-21E	2	15	576×419 / 22.7×16.5	EDA8137
Incubators	IB-05G	2	8	366×316 / 14.4×12.4	EDA8136
	IB-15G	2	12	446×366 / 17.6×14.4	EDA8135
	IB-25G	2	12	576×419 / 22.7×16.5	EDA8137
	IB-01E	2	11	366×316 / 14.4×12.4	EDA8136
	IB-11E	2	14	466×490 / 17.6×19.3	EDA8220
	IB-21E	2	15	566×466 / 22.3×17.6	EDA8219
	IL-11/ 11A	2	12	566×466 / 22.3×17.6	EDA8219
	IL-21/ 21A	3	26	466×490 / 17.6×19.3	EDA8220
	ILP-02	1	6	195×310 / 7.7×12.2	EDA8223
ILP-12	2	11	305×299 / 12.0×11.8	EDA8221	

Perforated Shelves

Stainless steel shelf system is removed easily without using tools for easy cleaning.

Product group	Models	Description		Cat. No.
		Maximum	W×L (mm/inch)	
Ovens	OF-02G	8	366×316 / 14.4×12.4	AAA12531
	OF-12G	12	446×366 / 17.6×14.4	AAA12532
	OF-22G	14	516×416 / 20.3×16.4	AAA12533
	ON-02G	8	366×316 / 14.4×12.4	AAA12531
	ON-12G	12	446×366 / 17.6×14.4	AAA12532
	ON-22G	12	576×419 / 22.7×16.5	AAA12535
	OF/ ON-01E	10	365×310 / 14.4×12.4	AAA14501
	OF/ ON-11E	14	445×360 / 17.6×14.4	AAA14502
	OF/ ON-21E	15	575×410 / 22.7×16.5	AAA14503
Incubators	IB-05G	8	366×316 / 14.4×12.4	AAA12531
	IB-15G	12	446×366 / 17.6×14.4	AAA12532
	IB-25G	12	576×419 / 22.7×16.5	AAA12535
	IB-01E	11	366×316 / 14.4×12.4	AAA12531
	IB-11E	14	466×490 / 17.6×19.3	AAA22522
	IB-21E	15	566×466 / 22.3×17.6	AAA22521
	IL-11/ 11A	12	566×466 / 22.3×17.6	AAA22521
	IL-21/ 21A	26	466×490 / 17.6×19.3	AAA22522
	ILP-02	6	195×310 / 7.7×12.2	AAA24501
ILP-12	11	305×299 / 12.0×11.8	AAA24502	

Mehrgan
Parto
Shar



Perforated Shelves

Stacking Kits

Stackable 2 steps up to maximize space efficiency with exclusive stacking kits. (Only for IB-G, IB-E series).

Product group	Models	Description	Cat. No.
Incubators	IB-05G	Stacking kits for stackable chambers	AAA21521
	IB-15G	Stacking kits for stackable chambers	AAA21522
	IB-25G	Stacking kits for stackable chambers	AAA21523
	IB-01E	Stacking kits for stackable chambers	AAA26501
	IB-11E	Stacking kits for stackable chambers	AAA26502
	IB-21E	Stacking kits for stackable chambers	AAA26503



Stackable Chamber

Tempered Viewing Windows

Optimized sample monitoring with threefold, tempered safety glass window. (Only for OF-G, ON-G series).

Product group	Models	Description	Cat. No.
		Thickness, W×L (mm/inch)	
Ovens	OF/ ON-02G	5, 150×280 / 0.2, 5.9×11.0	AAA12611
	OF/ ON-12G, 22G	5, 150×380 / 0.2, 5.9×15.0	AAA12612



Tempered Viewing Windows

Fan-Speed Adjuster

Adjustable inner air circulation speed, which protects spreading out powder sample. (Only for OF-G series).

Product group	Models	Description	Cat. No.
Ovens	OF/02G, 12G, 22G	Analog dial control	AAA12701

Mehrgan.
Parto
Shar



Fan Speed Adjuster